

S.NO	TITLE	PAGE.NO
1	Project Demonstration	03
2	Project Documentation	03
3	Feedback and Final Adjustments	04
4	Final Project Report Submission	05
5	Project Handover and Future Works	06

Phase 5: Project Demonstration & Documentation

Title: AI-Enhanced ICU Monitoring and Predictive Care System

Abstract:

The AI-Enhanced ICU Monitoring and Predictive Care System aims to modernize intensive care monitoring through a browser-based platform that simulates real-time vitals, predicts critical risks, and supports clinician decision-making. Built with HTML, CSS, and JavaScript in its prototype phase, the system integrates smart alerting, visual dashboards, and training simulations. Designed for remote monitoring and scalability, it prepares for future AI model deployment (LSTM/GRU) and IoT data integration. This document covers project demonstration, source code insights, system architecture, testing, feedback, and future works.

Index

1. Project Demonstration
2. Project Documentation
3. Feedback and Final Adjustments
4. Final Project Report Submission

5. Project Handover and Future Works
6. Screenshots (code + output)

1. Project Demonstration

Overview:

This section highlights the live demonstration of the ICU monitoring system. The purpose of this demonstration was to validate the core functionalities of the system, including its real-time monitoring capabilities, user interface responsiveness, intelligent alert logic, and training simulation module. It serves as a proof-of-concept for how a fully digital ICU environment can operate without physical sensors or hospital-grade hardware.

Demonstration Details:

- **System Walkthrough:** Navigation through dashboards showing heart rate, SpO₂, and warning alerts.
- **AI Risk Simulation:** Threshold-based logic simulates AI decision-making for vitals.
- **Training Engine:** Embedded ICU simulation video and dummy patient profiles used in the training module.
- **Performance Metrics:** Load tested with multiple patient cards, fast alert triggers under 2s latency.
- **Security Design:** Placeholder JWT structure and HTTPS preparation shown in the architecture.

Outcome:

The demonstration clearly showed that the ICU dashboard can function as an effective real-time monitoring and training tool. Stakeholders appreciated the intuitive design, quick alert generation, and the educational value of the training module. The system proved capable of simulating real ICU operations in a browser environment, confirming its readiness for backend integration, AI deployment, and institutional use in the next phase.

2. Project Documentation

Overview:

This section encompasses the technical foundation and complete blueprint of the ICU Monitoring System. The objective of this documentation is to provide clarity to developers, clinicians, IT administrators, and stakeholders on the system's internal structure, capabilities, and interaction flow. Every function, layout decision, and simulated logic component is documented for maintainability and future scalability. Documentation is critical for project continuity, handover, compliance, and iterative improvements.

Documentation Sections:

- **System Architecture:** HTML/CSS/JS dashboard with simulated data generators and alert triggers.
- **Code Documentation:** Source code files (e.g., [index.html](#), [style.css](#), [script.js](#)) structured and commented.
- **User Guide:** End-user views patient vitals, receives alerts; training page offers video-based learning.
- **Admin Guide:** Placeholder logic included for login roles; expansion planned in Phase 6.
- **Testing Reports:** Browser-based performance tests, UI responsiveness, alert logic validation.

Outcome:

The project is thoroughly documented, ensuring any developer or medical trainer can understand, modify, or enhance the system. It supports future transitions to AI integration or real patient monitoring.

3. Feedback and Final Adjustments

Overview:

This stage focuses on real-world usability testing and incorporating end-user feedback into the system. It emphasizes improving accessibility, interface clarity, and cross-device support based on evaluations by domain experts and test users. Feedback collection is a crucial part of software quality assurance and helps ensure the system meets practical healthcare needs and user expectations. This phase also validates improvements through final round testing, providing confidence for broader deployment.

Steps Taken:

- **Feedback Collection:** Received input from 2 test users and 1 mentor on alert clarity and screen layout.
- **Refinements:** Adjusted font sizes, color contrast, and alert positioning based on feedback.
- **Final Testing:** Verified all fixes across multiple browsers and patient profiles.

Outcome:

The system interface was enhanced based on user suggestions, improving color contrast, alert visibility, and mobile responsiveness. Test results confirmed smoother navigation and clearer user experience across browsers and screen sizes. These refinements increased confidence in the system's usability and user satisfaction.

4. Final Project Report Submission

Overview:

This section consolidates all efforts across Phases 1 through 5 into a structured report. The purpose of this report is to communicate the full lifecycle of the project—right from identifying ICU challenges, proposing AI-driven innovations, building and testing a prototype, and finally demonstrating its effectiveness. This documentation is essential for academic, institutional, and technical presentation of the project's credibility and impact.

Report Sections:

- **Executive Summary:** ICU dashboard simulates AI-based vitals tracking with future AI integration.
- **Phase Breakdown:** Covers design (P1), innovation (P2), implementation (P3), optimization (P4), and this final report (P5).
- **Challenges & Solutions:**
 - Alert fatigue: Solved via color-coded, prioritized alerts.

- Real-time simulation: Solved via interval-based JS and DOM optimizations.
- Data security planning: Role-based access and encryption placeholders included.

Outcomes:

The final report successfully documents the complete development lifecycle of the project. It clearly outlines the system's design, implementation, and demonstration, reinforcing its relevance in healthcare innovation. The report is ready for academic and institutional review.

5. Project Handover and Future Works

Overview:

This stage marks the transition from development to deployment preparation. It outlines the artifacts delivered to stakeholders, such as source code, system manuals, and training resources. It also defines a roadmap for enhancing and scaling the system. The goal is to ensure that future teams can expand the system to include real-time patient integration, full AI prediction, and wider institutional use. The overview acts as a strategic bridge between proof-of-concept and production-readiness.

Handover Details:

- **Next Steps:**

1. Add backend API with real-time vitals using Flask or Firebase.
2. Train and integrate LSTM model for risk scoring.
3. Add full authentication and role management.
4. Host on Vercel/Render with HTTPS.
5. Expand to mobile and multilingual support.

Outcome:

All project assets, including code, manuals, and documentation, are ready for handover. The system is now prepared for integration with backend APIs and AI models. Future plans include cloud deployment, authentication, and institutional trials.

SOURCE CODE:

HTML:

```
<!doctype html>
<html lang="en">
  <head>
    <meta charset="UTF-8" />
    <link rel="icon" type="image/svg+xml" href="/vite.svg" />
    <meta name="generator" content="Hostinger Horizons" />
    <meta name="viewport" content="width=device-width, initial-scale=1.0" />
    <title>Medical ICU ML Dashboard</title>
  </head>
  <body>
    <div id="root"></div>
    <script type="module" src="/src/main.jsx"></script>
  </body>
</html>
```

CSS:

```
@tailwind base;
@tailwind components;
@tailwind utilities;

/* === Base Theme Variables === */
@layer base {
  :root {
    --radius: 0.75rem;
    --background: 210 40% 98%;
    --foreground: 222.2 84% 4.9%;
    --card: 0 0% 100%;
    --card-foreground: 222.2 84% 4.9%;
    --popover: 0 0% 100%;
    --popover-foreground: 222.2 84% 4.9%;
    --input: 214.3 31.8% 91.4%;
    --ring: 262 83.3% 57.8%;
  }

  .dark {
    --background: 222.2 84% 4.9%;
    --foreground: 210 40% 98%;
    --card: 217.2 32.6% 17.5%;
    --card-foreground: 210 40% 98%;
    --popover: 222.2 84% 4.9%;
    --popover-foreground: 210 40% 98%;
  }
}
```

```

    --primary: 263.4 70% 50.4%;
    --accent: 191 91% 37%;

}

* {
    @apply border-border;
}

body {
    @apply bg-background text-foreground;
    font-feature-settings: "rlig" 1, "calt" 1;
    background: linear-gradient(135deg, hsl(var(--background)),
hsl(var(--muted)));
}

/* === Utilities === */
.glass-card {
    background: rgba(255, 255, 255, 0.1);
    backdrop-filter: blur(12px);
    border: 1px solid rgba(255, 255, 255, 0.2);
    box-shadow: 0 8px 32px rgba(31, 38, 135, 0.15);
}
.dark .glass-card {
    background: rgba(17, 25, 40, 0.75);
    border: 1px solid rgba(255, 255, 255, 0.125);
}

.gradient-bg {
    background: linear-gradient(135deg, hsl(var(--primary)), hsl(var(--accent)));
}

}

.card-hover {
    transition: all 0.4s cubic-bezier(0.4, 0, 0.2, 1);
}
.card-hover:hover {
    transform: translateY(-5px) scale(1.02);
    box-shadow: 0 20px 40px rgba(0, 0, 0, 0.15);
}

.pulse-animation {

```

```

    animation: pulse 2s cubic-bezier(0.4, 0, 0.6, 1) infinite;
  }
  @keyframes pulse {
    0%, 100% {
      opacity: 0.8;
      transform: scale(1);
    }
    50% {
      opacity: 1;
      transform: scale(1.05);
    }
  }
}

.grid-dashboard {
  display: grid;
  grid-template-columns: repeat(auto-fill, minmax(300px, 1fr));
  gap: 1.5rem;
}
@media (min-width: 1024px) {
  .grid-dashboard {
    grid-template-columns: repeat(auto-fill, minmax(350px, 1fr));
  }
}

.shimmer {
  background: linear-gradient(90deg, transparent 0%, rgba(255,255,255,0.2) 50%,
transparent 100%);
  background-size: 200% 100%;
  animation: shimmer 1.5s infinite;
}
@keyframes shimmer {
  0% { background-position: -200% 0; }
  100% { background-position: 200% 0; }
}

.input-gradient {
  background: linear-gradient(135deg, hsl(var(--background)), hsl(var(--muted)));
}
.input-gradient:focus {
  background: linear-gradient(135deg, hsl(var(--background)), hsl(var(--card)));
}

```

REACT:


```

import React, { useState, useMemo } from "react";
import {
  Card,
  CardContent,
  CardHeader,
  CardTitle,
} from "@components/ui/card";
import { Select, SelectTrigger, SelectValue, SelectContent, SelectItem } from
"@components/ui/select";
import { LineChart, Line, XAxis, YAxis, CartesianGrid, Tooltip, PieChart, Pie, Cell,
BarChart, Bar } from "recharts";
import { motion } from "framer-motion";
import { users } from "@data/users";

const getRiskColor = (risk: number) => {
  if (risk > 0.6) return 'bg-red-500';
  if (risk > 0.3) return 'bg-yellow-500';
  return 'bg-green-500';
};

const usePatientData = (selectedPatient: any) => {
  const avgRiskScore = useMemo(() => {
    const total = users.reduce((sum, u) => sum + u.sepsisRisk, 0);
    return (total / users.length).toFixed(2);
  }, []);

  const highRiskCount = useMemo(() => {
    return users.filter(u => u.sepsisRisk > 0.7).length;
  }, []);

  const conditionData = useMemo(() => {
    const conditionMap: Record<string, number> = {};
    users.forEach(u => {
      u.conditions.forEach(c => {
        conditionMap[c] = (conditionMap[c] || 0) + 1;
      });
    });
    return Object.entries(conditionMap).map(([name, count]) => ({ name, count }));
  }, []);

  const riskDistribution = useMemo(() => {
    return [

```

```

    { name: "Low", value: users.filter(u => u.sepsisRisk <= 0.3).length },
    { name: "Medium", value: users.filter(u => u.sepsisRisk > 0.3 && u.sepsisRisk
<= 0.6).length },
    { name: "High", value: users.filter(u => u.sepsisRisk > 0.6).length },
  ];
}, []);

return { avgRiskScore, highRiskCount, conditionData, riskDistribution };
};

```

```

const colors = ["#4ade80", "#facc15", "#f87171"];

```

```

export default function Analytics() {
  const [selectedPatientId, setSelectedPatientId] = useState("1");
  const selectedPatient = users.find(u => u.id === selectedPatientId)!;

  const { avgRiskScore, highRiskCount, conditionData, riskDistribution } =
    usePatientData(selectedPatient);

  return (
    <motion.div initial={{ opacity: 0 }} animate={{ opacity: 1 }} className="grid
grid-cols-1 gap-6 md:grid-cols-2 xl:grid-cols-3 p-4">
      <Card className="col-span-1">
        <CardHeader>
          <CardTitle>Average Risk Score</CardTitle>
        </CardHeader>
        <CardContent className="text-4xl font-bold">{avgRiskScore}</CardContent>
      </Card>

      <Card className="col-span-1">
        <CardHeader>
          <CardTitle>High Risk Patients</CardTitle>
        </CardHeader>
        <CardContent className="text-4xl font-bold">{highRiskCount}</CardContent>
      </Card>

      <Card className="col-span-1">
        <CardHeader>
          <CardTitle>Patient Selector</CardTitle>
        </CardHeader>
        <CardContent>
          <Select value={selectedPatientId} onChange={setSelectedPatientId}>
            <SelectTrigger className="w-full">

```

```

    <CardHeader>
    <CardTitle>Sepsis Risk Trend</CardTitle>
  </CardHeader>
  <CardContent>
    <LineChart width={300} height={200} data={selectedPatient.trendData}>
      <CartesianGrid strokeDasharray="3 3" />
      <XAxis dataKey="day" />
      <YAxis domain={[0, 1]} />
      <Tooltip />
      <Line type="monotone" dataKey="risk" stroke="var(--chart-primary)"
strokeWidth={2} dot={{ r: 4 }} />
    </LineChart>
  </CardContent>
</Card>

<Card className="col-span-1 md:col-span-2 xl:col-span-3">
  <CardHeader>
    <CardTitle>Patient Details</CardTitle>
  </CardHeader>
  <CardContent className="space-y-4">
    <div className="text-lg font-semibold">{selectedPatient.name}</div>
    <div className="text-sm text-muted-foreground">Age:
{selectedPatient.age}</div>
    <div className="text-sm text-muted-foreground">Gender:
{selectedPatient.gender}</div>
    <div className="text-sm text-muted-foreground">Conditions:
{selectedPatient.conditions.join(", ")}</div>
    <div className="text-sm text-muted-foreground">Sepsis Risk:
{selectedPatient.sepsisRisk}</div>
    <div className="h-2.5 rounded-full w-full">
      <div className={` ${getRiskColor(selectedPatient.sepsisRisk)} h-2.5
rounded-full`} style={{ width: `${selectedPatient.sepsisRisk * 100}%` }} />
    </div>
  </CardContent>
</Card>
</motion.div>
);
}

```

JAVA SCRIPT:

```

import path from 'node:path';
import react from '@vitejs/plugin-react';
import { defineConfig } from 'vite';

```

```

    )) {
      const shadow = node.shadowRoot;
      if (!shadow) return;
      const msg = shadow.querySelector('.message-body')?.textContent?.trim()
|| ";

      const file = shadow.querySelector('.file')?.textContent?.trim() || "";
      window.parent.postMessage({
        type: 'horizons-vite-error',
        error: msg + (file ? ' File: ' + file : "")
      }, '*');
    }
  }
});
observer.observe(document.documentElement, { childList: true, subtree: true
});
`,
injectTo: 'head'
},
{
  tag: 'script',
  attrs: { type: 'module' },
  children: `
    const origError = console.error;
    console.error = (...args) => {
      origError(...args);
      const err = args.find(a => a instanceof Error);
      const msg = err ? (err.stack || err.message) :
        args.map(a => typeof a === 'object' ? JSON.stringify(a) :
String(a)).join(' ');
      window.parent.postMessage({ type: 'horizons-console-error', error: msg },
'*');
    };
  `,
  injectTo: 'head'
},
{
  tag: 'script',
  attrs: { type: 'module' },
  children: `
    const origFetch = window.fetch;
    window.fetch = (...args) => {
      const url = args[0] instanceof Request ? args[0].url : args[0];

```

```

    if (url.startsWith('ws')) return origFetch(...args);
    return origFetch(...args)
      .then(async res => {
        if (!res.ok && !res.headers.get('content-type')?.includes('text/html')) {
          const text = await res.clone().text();
          console.error(`Fetch error from \${res.url}: \${text}`);
        }
        return res;
      })
      .catch(err => {
        if (!url.endsWith('.html')) console.error(err);
        throw err;
      });
  },
  injectTo: 'head'
}
];
return { html: "", tags: scripts };
}
}
];

```

```

// Simplified Vite config
export default defineConfig({
  plugins: [react(), ...injectScripts],
  server: {
    cors: true,
    headers: {
      'Cross-Origin-Embedder-Policy': 'credentialless',
    },
    allowedHosts: true,
  },
  resolve: {
    alias: {
      '@': path.resolve(__dirname, './src'),
    },
    extensions: ['.js', '.ts', '.jsx', '.tsx', '.json'],
  },
});

```

PYTHON:

```

project_management_tool = {
  "name": "project-management-tool",
  "version": "0.0.0",
  "scripts": {
    "dev": "vite",
    "build": "vite build",
    "preview": "vite preview"
  },
  "dependencies": {
    "@headlessui/react": "^1.7.20",
    "@heroicons/react": "^2.1.1",
    "@radix-ui/react-icons": "^1.3.0",
    "clsx": "^2.1.0",
    "lucide-react": "^0.358.0",
    "react": "^18.2.0",
    "react-dom": "^18.2.0",
    "react-icons": "^4.12.0",
    "shadcn/ui": "latest"
  },
  "devDependencies": {
    "@vitejs/plugin-react": "^4.2.1",
    "autoprefixer": "^10.4.17",
    "postcss": "^8.4.31",
    "tailwindcss": "^3.4.1",
    "vite": "^5.0.8"
  }
}

```

```

medical_icu_ml_dashboard = {
  "name": "medical-icu-ml-dashboard",
  "version": "0.0.0",
  "lockfileVersion": 3,
  "requires": true,
  "packages": {
    "": {
      "name": "medical-icu-ml-dashboard",
      "version": "0.0.0",
      "dependencies": {
        "@radix-ui/react-alert-dialog": "^1.0.5",
        "@radix-ui/react-avatar": "^1.0.3",
        "@radix-ui/react-checkbox": "^1.0.4",
        "@radix-ui/react-dialog": "^1.0.5",
        "@radix-ui/react-dropdown-menu": "^2.0.5",


```


```

"@radix-ui/react-label": "^2.0.2",
"@radix-ui/react-progress": "^1.0.3",
"@radix-ui/react-select": "^1.2.2",
"@radix-ui/react-separator": "^1.0.3",
"@radix-ui/react-slider": "^1.1.2",
"@radix-ui/react-slot": "^1.0.2",
"@radix-ui/react-tabs": "^1.0.4",
"@radix-ui/react-toast": "^1.1.5",
"@supabase/supabase-js": "^2.39.0",
"class-variance-authority": "^0.7.0",
"clsx": "^2.0.0",
"framer-motion": "^10.16.4",
"lucide-react": "^0.285.0"
}
}
}
}

```

Output Display Preview :





Medical ICU Dashboard

Sign in to access the dashboard

Username

Enter your username

Password

Enter your password

[Sign In](#)

Logged out
You have been successfully logged out.

Medical ICU ML Dashboard

Patients Management

Patient Search & Filter

...

Add New Patient

Enter the details for the new patient. Click save when you're done.

Name

arun

Age

74

Condition

Primary medical condition

Status

Stable

Vitals

Heart Rate (bpm)

e.g., 70

Blood Pressure (mmHg)

e.g., 120/80

Temperature (°C)

e.g., 37.0

SpO2 (%)

e.g., 98

Respiratory Rate (bpm)

e.g., 16


Cancel

Save Patient

Medical ICU ML Dashboard

Search patients, metrics, or ML models...

<div><div>JO</div><div>John Doe</div><div>ID: P101 Age: 45</div><div>Stable</div></div> <div><div>Condition: Pneumonia</div><div><div>HR: 75 bpm</div><div>BP: 120/80</div><div>Temp: 37°C</div><div>SpO2: 97%</div><div>RR: 18 bpm</div></div><div><div>Last Update: 1h ago</div><div>View</div></div></div>	<div><div>JA</div><div>Jane Smith</div><div>ID: P102 Age: 62</div><div>Critical</div></div> <div><div>Condition: Heart Failure</div><div><div>HR: 95 bpm</div><div>BP: 100/60</div><div>Temp: 36.5°C</div><div>SpO2: 92%</div><div>RR: 24 bpm</div></div><div><div>Last Update: 15m ago</div><div>View</div></div></div>	<div><div>RO</div><div>Robert Johnson</div><div>ID: P103 Age: 50</div><div>Improving</div></div> <div><div>Condition: Sepsis</div><div><div>HR: 88 bpm</div><div>BP: 110/70</div><div>Temp: 38.2°C</div><div>SpO2: 95%</div><div>RR: 20 bpm</div></div><div><div>Last Update: 30m ago</div><div>View</div></div></div>	<div><div>EM</div><div>Emily White</div><div>ID: P104 Age: 33</div><div>Stable</div></div> <div><div>Condition: Asthma Exacerbation</div><div><div>HR: 80 bpm</div><div>BP: 115/75</div><div>Temp: 36.8°C</div><div>SpO2: 99%</div><div>RR: 16 bpm</div></div><div><div>Last Update: 2h ago</div><div>View</div></div></div>
<div><div>AR</div><div>arun</div><div>ID: P391 Age: 74</div><div>Improving</div></div> <div><div>Condition: 23</div><div><div>HR: 56 bpm</div><div>BP: 53/44</div><div>Temp: 34°C</div><div>SpO2: 42%</div><div>RR: 22 bpm</div></div><div><div>Last Update: Just now</div><div>View</div></div></div>			



Medical ICU ML Dashboard

Search patients, metrics, or ML models...

A

All Types

All Statuses

Discharge Summary

Completed

For John Doe (ID: P101)

Date: 2025-05-10

Physician: Dr. Smith

Download

Progress Note

Pending Review

For Jane Smith (ID: P102)

Date: 2025-05-12

Physician: Dr. Emily Jones

Download

Sepsis Alert Report

Completed

For Robert Johnson (ID: P103)

Date: 2025-05-08

Physician: Dr. Brown

Download

Consultation Note

Generated

For Emily White (ID: P104)

Date: 2025-05-11

Physician: Dr. Davis

Download

Lab Results Summary

Completed

For John Doe (ID: P101)

Date: 2025-05-09

Physician: Lab System

Download

